

# AlbuRich P140 | PRP | P15

SKU 67675 (P140) | 67685 (PRP) | 67645 (P15)

## SCIENTIFIC ANALYSIS

ATTRIBUTE	SPECIFICATION	TYPICAL
Protein (dry basis)	≥96.0%	≥99.0%
Purity (albumin)	≥98%	100%
Moisture	≤5.0%	<2.0%
pH	6.5 to 7.5	6.8 to 7.2
BIGG	None Detected	None Detected
Lipid Content	≥ 1–3 µg Total/mg	≥ 2.5 µg Total/mg
Protease	<0.001 units/mg	<0.001 units/mg
Microbial Content	≤100 cfu/g	≤100 cfu/g
Endotoxin	≤3 EU/mg	≤0.1 EU/mg

## PRODUCT BENEFITS

- High Purity, Low Interference, Low Background
- No Detectable IgG or Protease Activity
  - Ensures assay integrity, component stability, and cell culture suitability
- Superior Solubility/Filterability for Ease of Use

## PROCESS BENEFITS

- Origin Traceability to Facilitate Regulatory Approval
- Closed-Loop Manufacturing\*
  - Eliminates the possibility of contamination and maximizes reproducibility
- Product Development and Technical Services designed for maximum responsiveness and flexibility

## QUALIFICATIONS

- Manufactured from MPI-inspected bovine plasma
- Validated for Prion (TSE) Clearance
- USDA/APHIS Registered for Shipment to the EU
- Manufactured according to cGMPs
- TSE Certificate of Suitability

ONLY UNDER 4C AND CHAPTER 20

Continued on the other side.



2425 SE Oak Tree Court • Ankeny, Iowa 50021 • USA  
phone 515.289.5170  
phb1.com

## DESCRIPTION

Proliant Bovine Serum Albumin (BSA) is a high purity, lyophilized “Fraction V” powder derived from the plasma of MPI inspected healthy animals from abattoirs located in New Zealand.

AlbuRich was developed with the understanding that there are a number of proprietary media formulations in the cell culture market. The AlbuRich formulations were created to help you find what works best with your current media formulations and are available in a variety of options depending on fatty acid required for your applications.

## PHYSICAL CHARACTERISTICS

Proliant BSA is a highly soluble white-to-yellow with tan-to-green cast lyophilized powder. A solution of up to 30% in de-ionized water is clear-to-slightly hazy and virtually particulate-free.

## PACKAGING

Standard package sizes are 10 kg and 1 kg. Sample sizes of 100 g are available for evaluation.

PACKAGE	SIZE
Standard	10 kg
	1 kg
Sample	100 g

## STORAGE RECOMMENDATIONS

Store at 2 to 8°C for up to 3 years from the date of manufacture

## ORDERING INFORMATION

Please order by referencing:

10 kg:	67675 (P140)	1 kg:	67675-1kg (P140)
	67685 (PRP)		67685-1kg (PRP)
	67645 (P15)		67645-1kg (P15)

\*Proliant's unique **Closed-Loop** Manufacturing means that we collect only in exclusive abattoirs with our own USDA/MPI approved proprietary collection system. No other manufacturer can claim this important closed link, as they must buy their raw material from outside sources and the spot market. What's more, from the time our product enters our tankers to the moment it is loaded into the dryers, the product is not exposed to the external environment. Once the raw material reaches our cGMP plant, the material is loaded directly into our process and manufactured through a series of closed tanks, lines, filters and separation equipment. In our plant there are no exposed filter presses or open tanks, nor is any non-food grade or noxious solvent material introduced into the process. There is also no possibility of contamination with other animal-derived products or bovine tissues in our completely dedicated facility.