

Bovine Serum Albumin

Caprylate Stabilized (Standard Grade) pH 5.2

SKU 68500

SCIENTIFIC ANALYSIS

ATTRIBUTE	SPECIFICATION	TYPICAL
Protein (dry basis)	≥98.0%	≥99.0%
Purity (albumin)	≥98%	100%
Solubility (4% water)	Clear-slightly hazy	Clear-slightly hazy
Moisture	≤5.0%	<2.0%
pH (10% water)	5.0 to 5.6	5.0 to 5.4
Ash	<2%	<2%
Protease	<0.001 unit/g	≤0.005 Units/mg
Endotoxin	≤1 EU/mg	≤0.1 EU/mg
Heavy Metals (Pb)	≤10 ppm	<1 ppm

Na, Cl, K, P, Ca are tested and can be reported on the Certificate of Analysis if requested.

PRODUCT BENEFITS

- High Purity, Low Interference, Low Background
- No Detectable IgG, Endotoxin, Enzyme or Protease Activity
 - Ensures assay integrity, component stability, and cell culture suitability
- Superior Solubility/Filterability for Ease of Use

PROCESS BENEFITS

- Origin Traceability to Facilitate Regulatory Approval
- Closed-Loop Manufacturing*
 - Eliminates the possibility of contamination and maximizes reproducibility
- Product Development and Technical Services designed for maximum responsiveness and flexibility

QUALIFICATIONS

- Manufactured from USDA-inspected bovine plasma
- Validated for Prion (TSE) Clearance
- USDA/APHIS Registered for Shipment to the EU
- Manufactured according to FDA cGMPs
- TSE Certificate of Suitability

Continued on the other side.



2425 SE Oak Tree Court • Ankeny, Iowa 50021 • USA
phone 515.289.5170
phb1.com

DESCRIPTION

Proliant Bovine Serum Albumin (BSA) is a high purity, lyophilized “Fraction V” powder derived from the plasma of USDA/MPI inspected healthy animals from abattoirs located in the United States and New Zealand.

Proliant Caprylate Stabilized (Standard Grade) pH 5.2 BSA is the established standard for lower pH applications because of its undetectable protease and IgG along with the fact that it is also free of contaminants that can cause background interference. At pH 5.2, our BSA has optimal antifouling properties, making it a superior blocking agent across a variety of applications.

PHYSICAL CHARACTERISTICS

Proliant BSA is a highly soluble white-to-yellow with tan-to-green cast lyophilized powder. A solution of up to 30% in de-ionized water is clear-to-slightly hazy and virtually particulate-free.

PACKAGING

Standard package sizes are 10 kg and 1 kg. Sample sizes of 100 g are available for evaluation.

PACKAGE	SIZE
Standard	10 kg
	1 kg
Sample	100 g

STORAGE RECOMMENDATIONS

Store sealed in cool, dry environment for 5 years.

ORDERING INFORMATION

Please order by referencing:

10 kg:	68500
1 kg:	68500-1kg

*Proliant's unique **Closed-Loop** Manufacturing means that we collect only in exclusive abattoirs with our own USDA approved proprietary collection system. No other manufacturer can claim this important closed link, as they must buy their raw material from outside sources and the spot market. What's more, from the time our product enters our tankers to the moment it is loaded into the dryers, the product is not exposed to the external environment. Once the raw material reaches our cGMP plant, the material is loaded directly into our process and manufactured through a series of closed tanks, lines, filters and separation equipment. In our plant there are no exposed filter presses or open tanks, nor is any non-food grade or noxious solvent material introduced into the process. There is also no possibility of contamination with other animal-derived products or bovine tissues in our completely dedicated facility.